

Evco Controller Operating Manual

Thank you for your purchase. This manual contains: set-up instructions, unit specs, safety information, controller operation, and maintenance steps.

K212, K225, K230, K249, K272



IMPORTANT:

Your unit is preprogrammed

Place your unit in the desired location. Plug the unit in and allow it to cool and become stable for a minimum of 24 hours before logging temperature or stocking products.

Be careful when setting or changing temperatures

WARNING: Changing some controller parameters can damage your unit and/or result in a loss of product. K2 will not be held responsible for losses due to unauthorized parameter changes.

Changing advanced parameters may damage the unit or void your warranty. Please contact K2 before attempting to change advanced parameters.

Unfamiliar with the operation of a K2 controller?

Use the video tutorials on our website resources page or call us for assistance with special parameters.





CONTENTS

- 1. Front Panel Commands
- 2. General Operations
- 3. Parameters

- 4. Alarms
- 5. Electrical Connection
- 6. Temperature Logging

1 Front Panel Commands





LED	MODE	MEANING		
☀	-ON -OFF -Blinking	-Compressor ON -Compressor OFF -Compressor protection activated/set point temperature menu		
脊	-ON - Active Defrost/ pre-dripping cycle -OFF -No action. -Blinking -Defrost delay time/active dripping cycle			
@	-ON -OFF -Blinking	-Evaporator fan ON -Evaporator fan OFF -Evaporator fan stop		
НАССР	-ON -OFF -Blinking	-HACCP alarm recorded in EVlink module -no action -no action		
٢	-ON -OFF -Blinking	-Active Energy saving mode -no action -no action		
عر	-ON -OFF -Blinking	-request of compressor's servicing -no action -active settings mode/ active access to additional functions/ active connection with EVlink module		
°C/°F	-ON -OFF -Blinking	-Normal temperature view -no action -Active overheating/overcooling cycle		
AUX	-ON -OFF -Blinking	-Auxiliary load ON -Auxiliary load OFF -Auxiliary Load activated by digital input / Auxiliary Load activation delay		
Ú	-ON -OFF -Blinking	-Device OFF -Device ON -Device ON/OFF mode		

2 General Operations

Keypad Unlocking

- The keypad will lock automatically after 30 seconds if the default parameter Loc=1 is set.
- To unlock the keypad, touch the SET key for 1 second. Then the display will show the label "**UnL**".

Temperature Settings

- If the keypad is locked, first unlock it.
- Touch the SET key then set the desired temperature by pressing the UP or DOWN keys within 15 seconds.
- Your set temperatures must be within temperature range high and low limits. The limits range of the set point (parameters **r1** and **r2**)
- Press SET to confirm.

Cabinet Light

- If the parameter is u1=0
- Touch the ON/STAND-BY key.

Buzzer

• If the parameters u1=3 and u4=1 touch any key to shut down the buzzer alarm.

Manual Defrost

- Unlock the keypad
- Press the UP key, holding it for 2 seconds.
- If the parameter **P3=1** and the evaporator temperature value is lower than the parameter **d2**, defrost cycle will start.

Compressor Opperational Time

- Make sure the keypad is unlocked and press the DOWN key for 4 seconds.
- Scroll through the menu labels with the UP and DOWN keys.
 - CH label: displaying compressor operating hours.
 - rCH label: compressor operating hours reset.
- **nS1** label: compressor star-up time.
- To access the label press SET.
- In order to reset the compressor operating hours once selected the **rCH** label, insert the password "**149**" using the UP and DOWN keys then confirm touching the SET key.
- •Touch the ON/STANDBY key to exit the procedure.

Displaying Temperature Probes

- Ensure the keypad is unlocked then press the DOWN key for 4 seconds.
- Scroll through the menu's labels by the UP or DOWN key.
- •Pb1: cabinet temperature probe

(if parameter **P4=0,1** or **2**); inlet air

- temperature probe (if parameter **P4=3**).
- •Pb2: Evaporator temperature probe
- (if parameter P3=1 or 2)
- •**Pb3:** Auxiliary temperature probe (If **P4=1, 2** or **3**).
- To access the label press SET.
- Touch the ON/STANDBY key to exit the procedure.



Alarm code	Code description	Solution
Pr1	Cabinet probe alarm	-Check the parameter P0.
Pr2	Evaporator probe alarm	-Check the status of the probe.
Pr3	Auxiliary probe alarm	-Check the electrical connection.
		-Replace the probe.
rtc	Date and time alarm	Set date, time and day of the week.
AL	Low Temperature Alarm	Check the parameters AA, A1 and A2
AH	High Temperature Alarm	Check the parameter AA, A4 and A5
id	Door open alarm	Check the parameter i0 and i1
PF	Power failure alarm	-Check electrical connection
		-Touch any key to shut the buzzer off
сон	High condenser warning	-Check if the condenser probe is
		installed.
		-Check the parameter C6
		-Check the condenser coil is clean.
CSd	High condensation alarm	Check if the condenser probe is
		installed.
		-Check the parameter C7.
		-Check the condenser coil is clean.
		-Reboot the device.
iA	Multi-function input alarm	-Check the parameters i5 and i6
Cth	Compressor thermal switch alarm	-Check the parameters i5 and i6
th	Global thermal switch alarm	-Check the parameters i5 and i6
		-Reboot the device
dFd	Defrost time out alarm	-Check the parameters d2, d3 and d11

4 Alarms

5 Electrical Connection





3 Parameters

WARNING: Do not modify parameters in this unit without fully understanding their function. Changing advanced parameters may damage the unit and void your warranty. Please contact K2 before attempting to change advanced parameters.

Accessing Operational Parameters

- Touch the SET key for 4 seconds, the monitor will display the label "PA"
- Press the SET key and insert the password "-19"
- Press SET key to confirm
- Scroll through the parameters list using the UP or DOWN key
- For modifying a parameter value, press SET key at the parameter label then adjust the value by the UP or DOWN key
- Press SET key to confirm the change
- Press SET key for 4 seconds or do not operate for 60 sec to exit the procedure.

	Ν	Parameter	Default	Setpoint	
₽≣	1	SP	0.0	setpoint	
0	8	P3	1	evaporator probe function	
÷.	13	r0	2.0	setpoint differential	
_	14	r1	-50	minimum setpoint	
	15	r2	50.0	maximum setpoint	
•	33	d0	8	automatic defrost interval	
	34	d1	0	defrost type	
	35	d2	8.0	threshold for defrost end	
	36	d3	30	defrost duration	
Å	51	AA	0	select value for high/low temperature alarms	
	52	A1	-10.0	threshold for low temperature alarm	
	53	A2	2	low temperature alarm type	
	54	A4	10.0	threshold for high temperature alarm	
	55	A5	2	high temperature alarm type	
ą	74	iO	5	door switch input section	
	75	i1	0	door switch input activation	
	76	i2	30	open door alarm delay	

Basic parameters

Full list of parameters

	Ν.	PAR.	DEF.	SETPOINT	MIN MAX.
9 -	1	SP	-20	setpoint	r1 r2
	Ν.	PAR.	DEF.	ANALOGUE INPUTS	MIN MAX.
	2	CA1	0.0	cabinet probe offset	-25 25 °C
	3	CA2	0.0	evaporator probe offset	-25 25 °C
	4	CA3	0.0	condenser probe offset	-25 25 °C
[5	PO	1	probe type	0 = PTC 1 = NTC
[6	P1	0	enable °C decimal point	0 = no 1 = yes
	8	P3	1	evaporator probe function	0 = disabled
O.					1 = defrost + fan
					2 = fan
	9	P4	2	configurable input function	0 = door switch input
					1 = condenser probe
					2 = condenser probe + door
	10	DE		value diaplayed	switch input
	10	42		value displayed	0 = cabinet temperature
					2 = evaporator temperature
					3 = condenser temperature
ł	11	P8	0	display refresh time	0 250 s : 10
I					
	N.	PAR.	DEF.	REGULATION	MIN MAX.
	N. 12	PAR. r0	DEF. 3.0	REGULATION setpoint differential	MIN MAX. 1 15 °C
	N. 12 13	PAR. r0 r1	DEF. 3.0 -30	REGULATION setpoint differential minimum setpoint	MIN MAX. 1 15 °C -99 °C r2
	N. 12 13 14	PAR. r0 r1 r2	DEF. 3.0 -30 -10	REGULATION setpoint differential minimum setpoint maximum setpoint	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C
J.	N. 12 13 14 15	PAR. r0 r1 r2 r4	DEF. 3.0 -30 -10 0.0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C
*	N. 12 13 14 15 16	PAR. r0 r1 r2 r4 r5	DEF. 3.0 -30 -10 0.0 0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling
*	N. 12 13 14 15 16	PAR. r0 r1 r2 r4 r5	DEF. 3.0 -30 -10 0.0 0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating
*	N. 12 13 14 15 16 17	PAR. r0 r1 r2 r4 r5 r6	DEF. 3.0 -30 -10 0.0 0 0.0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint setpoint	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C
*	N. 12 13 14 15 16 17	PAR. r0 r1 r2 r4 r5 r6	DEF. 3.0 -30 -10 0.0 0 0.0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C
*	N. 12 13 14 15 16 17 17	PAR. r0 r1 r2 r4 r5 r6 r7	DEF. 3.0 -30 -10 0.0 0 0.0 0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min
*	N. 12 13 14 15 16 17 17 18 19	PAR. r0 r1 r2 r4 r5 r6 r7 r8	DEF. 3.0 -30 -10 0.0 0 0.0 0 0 0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration DOWN key additional function	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min 0 = disabled
*	N. 12 13 14 15 16 17 17 18 19	PAR. r0 r1 r2 r4 r5 r6 r7 r8	DEF. 3.0 -30 -10 0.0 0 0.0 0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration DOWN key additional function	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min 0 = disabled 1 = overcooling/overheating 2 = opergy caprime
*	N. 12 13 14 15 16 17 18 19	PAR. r0 r1 r2 r4 r5 r6 r7 r8	DEF. 3.0 -30 -10 0.0 0 0.0 0 0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration DOWN key additional function	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min 0 = disabled 1 = overcooling/overheating 2 = energy saving
*	N. 12 13 14 15 16 17 18 19 20	PAR. r0 r1 r2 r4 r5 r6 r7 r8 r13	DEF. 3.0 -30 -10 0.0 0 0.0 0 0 25.0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration DOWN key additional function proportional band (relative to setpoint)	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min 0 = disabled 1 = overcooling/overheating 2 = energy saving 0 99 °C setpoint + r13
*	N. 12 13 14 15 16 17 18 19 20 21	PAR. r0 r1 r2 r4 r5 r6 r7 r8 r13 r14	DEF. 3.0 -30 -10 0.0 0 0.0 0 25.0	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration DOWN key additional function proportional band (relative to setpoint) integral action time	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min 0 = disabled 1 = overcooling/overheating 2 = energy saving 0 99 °C setpoint + r13 0 99 min
*	N. 12 13 14 15 16 17 18 19 20 21	PAR. r0 r1 r2 r4 r5 r6 r7 r8 r13 r14	DEF. 3.0 -30 -10 0.0 0 0.0 0 25.0 10	REGULATION setpoint differential minimum setpoint maximum setpoint setpoint offset in energy saving cooling or heating operation setpoint offset in overcool- ing/overheating overcooling/overheating duration DOWN key additional function proportional band (relative to setpoint) integral action time	MIN MAX. 1 15 °C -99 °C r2 r1 199 °C 0 99 °C 0 = cooling 1 = heating 0 99 °C 0 240 min 0 = disabled 1 = overcooling/overheating 2 = energy saving 0 99 °C setpoint + r13 0 99 min

	Ν.	PAR.	DEF.	COMPRESSOR	MIN MAX.
	23	CP0	0	time compressor at 85 Hz after	0 100 s x 10
				power-on	
	24	C0	1	compressor on delay after pow-	0 240 min
				er-on	
	25	C2	3	compressor off minimum time	0 240 min
	26	C3	0	compressor on minimum time	0 240 s
				(minimum speed)	
	27	C4	5	compressor off time during cabi-	0 240 min
				net probe alarm	
	28	C5	10	compressor on time (maximum	0 240 min
e				speed) during cabinet probe	
				alarm	
	29	C6	55.0	threshold for high condensation	0 199 °C
				warning	differential = 2 °C/4
	30	C7	60.0	threshold for high condensation	0 199 °C
				alarm	
	31	C8	1	high condensation alarm delay	0 15 min
	32	C9	5	consecutive time cabinet tem-	0 99 h
				perature in proportional band for	0 = disabled
				compressor at maximum speed	until cabinet temperature <
					setpoint
	33	C10	0	compressor hours for service	0 999 h x 10
		040	DEE		
	N.	PAR.	DEF.	DEFROST (IF r5 = 0)	MIN MAX.
	34	au	12	automatic derrost interval	
					d = 3 maximum interval
	35	d1	1	defrost type	0 = electric
			-	demost type	1 = hot gas
					2 = compressor stopped
	36	d2	6.0	threshold for defrost end	-99 99 °C
	37	d3	30	defrost duration	0 99 min
					se P3 = 1, maximum duration
	38	d4	0	enable defrost at power-on	0 = no 1 = yes
	39	d5	0	defrost dealy after power-on	0 99 min
	40	d6	1	value displayed during defrost	0 = cabinet temperature
					1 = display locked
					2 = dEF label
	41	d7	3	dripping time	0 15 min
	42	d8	0	defrost interval counting mode	0 = device on hours
					1 = compressor on hours
					2 = hours evaporator tem-
					perature < d9
					3 = adaptive
	43	d9	0.0	evaporation threshold for auto-	-99 99 °C
•				matic defrost interval counting	0
•	44	d11	1	enable derrost timeout alarm	0 = no 1 = yes
	45	a15	0	compressor on consecutive time	-20 99 min
-				for not gas defrost	n negative values, dripping
	46	d16	0	pre-dripping time for bot ass do	
	40	010		frost	0 99 11111
				11030	I

	46	d16	0	pre-dripping time for hot gas de-	0 99 min
				frost	
	47	d18	0	adaptive defrost interval	0 999 min
					if compressor on + evapora-
					tor temperature < d22
					0 = only manual
	48	d19	0.0	threshold for adaptive defrost	0 40 °C
				(relative to optimal evaporation	optimal evaporation tempera-
				temperature)	ture - d19
	49	d20	0	compressor on consecutive time	0 999 min
				for defrost	0 = disabled
	50	d21	0	compressor on consecutive time	0 500 min
				for defrost after power-on and	if (regulation temperature -
				overcooling	setpoint) > 10°C/20
					0 = disabled
	51	d22	0.0	evaporation threshold for adap-	-10 10 °C/°F
				tive defrost interval counting	optimal evaporation tempera-
				(relative to optimal evaporation	ture + d22
				temperature)	
	Ν.	PAR.	DEF.	ALARMS	MIN MAX.
	52	A1	0.0	threshold for low temperature	0 99 °C
				alarm (relative to setpoint)	0 = disabled
					cabinet temperature - A1
	53	A4	50.0	threshold for high temperature	0 99 °C
				alarm (relative to setpoint)	0 = disabled
					cabinet temperature + A4
	54	A6	12	high temperature alarm delay af-	0 99 min x 10
\mathbf{N}				ter power-on	
	55	A7	15	high/low temperature alarms de-	0 240 min
				lay	
	56	A8	15	high temperature alarm delay af-	0 240 min
				ter defrost	
	57	A9	15	high temperature alarm delay af-	0 240 min
				ter door closing	
	58	A11	2.0	high/low temperature alarms re-	1 15 °C
	N.	PAR.	DEF.	FANS	MIN MAX.
	59	FO	1	evaporator fan mode during	0 = off 1 = on
			-	normal operation	2 = according to F15 and
					F16 if compressor off, on
					if compressor on
					3 = thermore gulated (with
					F1)
					4 = thermoregulated (with
					F1) if compressor on

K SCIENTIFIC	
---------------------	--

	60	F1	0.1	threshold for evaporator fan op- eration	1 15 °C
	61	F2	0	evaporator fan mode during de- frost and dripping	$0 = off \qquad 1 = on$ 2 = according to F0
-	62	F3	2	evaporator fan stop maximum duration	0 15 min
S	63	F4	30	evaporator fan off time during energy saving	0 240 s x 10
	64	F5	30	evaporator fan on time during energy saving	0 240 s x 10
	65	F6	30	evaporator fan on time after compressor on	0 240 s x 10 if F0 = 3 or 4
	66	F7	20.0	threshold for evaporator fan on after dripping (relative to setpoint)	-99 99 °C setpoint + F7
	67	F8	2.0	threshold for evaporator fan op- eration differential	1 15 °C
	68	F9	10	evaporator fan off delay after compressor off	0 240 s if F0 = 2
	69	F10	0	evaporator fan and condenser fan off minimum time	0 240 s
	70	F11	10.0	threshold for condenser fan on	0 99 °C
	71	F12	0	condenser fan off delay after compressor off	0 240 s if P4 = 0
	72	F13	2.0	threshold for condenser fan on differential	1 15 °C
	73	F14	0	condenser fan mode	0 = thermoregulated (with F11) 1 = thermoregulated (with F11) if compressor on
	74	F15	60	evaporator fan off time with compressor off	0 240 s if F0 = 2
	75	F16	10	evaporator fan on time with compressor off	0 240 s if F0 = 2
	Ν.	PAR.	DEF.	DIGITAL INPUTS	MIN MAX.
	76	iO	2	door switch input function	 0 = disabled 1 = compressor + evaporator fan off 2 = evaporator fan off 3 = cabinet light on 4 = compressor + evaporator fan off, cabinet light on 5 = evaporator fan off + cabinet light on
	77	11	1	door switch input activation	0 = with contact closed 1 = with contact open
-	78	i2	0	open door alarm delay	-1 120 min -1 = disabled
	79	i3	-1	regulation inhibition maximum time with door open	-1 120 min -1 = until the closing
	80	110	0	door closed consecutive time for energy saving	0 999 min after regulation temperature < SP 0 = disabled

	81	i13	0	number of door openings for de-	0 240
				frost	0 = disabled
	82	i14	0	door open consecutive time for	0 240 min
				defrost	0 = disabled
	Ν.	PAR.	DEF.	DIGITAL OUTPUTS	MIN MAX.
	83	uc	1	enable relay K1 and relay K4 in- version	0 = no 1 = yes
	84	uc2	0	relay K2 configuration	0 = evaporator fan 1 = dripping heaters
	85	uc3	1	relay K3 configuration	0 = condenser fan 1 = cabinet light 2 = demisting 3 = on/stand-by 4 = compressor
*	86	u1	6	relay K4 configuration	0 = cabinet light 1 = demisting 2 = button-operated load 3 = alarm 4 = door heaters 5 = heater for neutral zone 6 = condenser fan 7 = on/stand-by
	87	u2	1	enable cabinet light and button- operated load in stand-by	0 = no 1 = yes manual
	88	u4	0	enable alarm output off silencing the buzzer	0 = no 1 = yes
	89	u5	-1.0	threshold for door heaters on	-99 99 °C differential = 2 °C/4
	90	u6	5	demisting on duration	1 100 min x 10
	91	u7	-5.0	neutral zone threshold for heat- ing (relative to setpoint)	-99 99 °C/°F differential = 2 °C/4 setpoint + u7
	Ν.	PAR.	DEF.	ENERGY SAVING (if $r5 = 0$)	MIN MAX.
<u>.</u>	92	HE2	0	energy saving maximum duration	0 999 min -1 = until the door opening
	93	HE3	0	consecutive time without operat- ing on keys for low consumption	0 240 min
~	N.	PAR.	DEF.	SAFETIES	MIN MAX.
$\overline{\mathbf{A}}$	94	POF	1	enable ON/STAND-BY key	0 = no 1 = yes
\sim	95	PAS	22	password	-99 999

6 Temperature Logging

SCIENTIFIC

When storing vaccines you may be required to preform a field validation test. A NIST callibrated external data logger can be used for this purpose. A data logger with text, email or online access is an added layer of protection for your product load in the event of a temperature excursion. K2 offers NIST calibrated data loggers to match your unit.

Service

K2 Scientific want to make sure you are happy with your purchase. There are several ways for you to contact us with questions or service needs. Be sure to include your four digit order number or at least your model number handy to speed up the process.

- 1; Contact us via our chat feature at www.k2sci.com
- 2: Email support@k2sci.com
- 3: Call 800-218-7613